Mode d'emploi  
Gebrauchsanweisung  
Instrucciones de utilización  
Návod k použití  
Gebruiksaanwijzing  
Bruksanvisning

Instructions for Use  
Istruzioni per l'uso  
Instruções de utilização  
Οδηγίες χρήσης  
Instrukcja użycia  
Kullanım kılavuzu

## DUTAU-NOVATECH®

Bronchoscope rigide Tuhý bronchoskop

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## 1 About this Document

### 1.1 Symbols Glossary

Symbol	Description
	Caution: Consult Instructions for Use
	Do not use if package is damaged
	Keep away from direct sunlight
	Keep dry
	Medical device
	Catalog number
	Batch code
	Unique Device Identification (UDI)
	Quantity per packaging unit
	Manufacturer
	Date of manufacture
	Consult Instructions for Use. The Instructions for Use are provided in electronic form (e-labelling).

Table 1: Symbols Glossary

### 1.2 Abbreviations

- WD: Washer Disinfector

### 1.3 Terminology

- Distal: From the surgeon's view further away
- Proximal: From the surgeon's view closer

### 1.4 Safety Information Marking

#### **WARNING**

Non-compliance may result in serious injuries, serious deterioration of the general condition or the death of the patient, user, or a third party.

#### **NOTICE**

Product damage or other damage may occur in case of non-compliance.

## 2 Important Safety Information

### **⚠ WARNING**

- Before using the product, read the Instructions for Use. Adhere to and save the Instructions for Use. Otherwise there are risks to the health of your patient.
- Do not modify the product. Otherwise there are risks to the health of your patient.

ATTENTION: In case that any serious incident has occurred in relation to the device the incident should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

## 3 Product Codes / REF

REF	Product
02BRD SET	DUTAU NOVATECH Rigid Bronchoscope SET
02BRD SET FR	DUTAU NOVATECH Rigid Bronchoscope SET
02BRD1100	DUTAU NOVATECH Bronchoscope head

Table 2: Product Codes

## 4 Scope of Delivery

### 4.1 Bronchoscope Head (REF 02BRD1100)

- 1 x Bronchoscope head corpus
  - 1 x ventilation element (including 2 x O-ring ID 16 x 1.25)
  - 1 x axial retainer (including 1 x O-ring ID 15 x 1.5)
  - 1 x axial cap
- 1 x Instructions for Use

### 4.2 Rigid Bronchoscope Set (02BRD SET / 02BRD SET FR)

REF	Product	Color	Ø inside / outside [mm]	Length [cm]
02BRD1100	Bronchoscope head			
02BRD1110 <sup>a</sup>	Injection cannula			
02BRD1120 <sup>a</sup>	Lateral double port (including 2 x lateral cap)			
NO117 <sup>a</sup>	Instructions for Use			
02BRD2514 <sup>a</sup>	Tracheal Tube, size 14	Blue and red	13 / 14	25

REF	Product	Color	Ø inside / outside [mm]	Length [cm]
02BRD3514 <sup>a</sup>	Bronchial Tube, size 14	Red	13 / 14	35
02BRD3512 <sup>a</sup>	Bronchial Tube, size 12	Black	10,7 / 12	35
02BRD3511 <sup>a</sup>	Bronchial Tube, size 11	Green	10,5 / 11,2	35
02BRD3509 <sup>a, b</sup>	Bronchial Tube, size 9	White	8,5 / 9,2	35
<sup>a</sup> Available separately as spare parts. <sup>b</sup> Not with 02BRD SET FR				

Table 3: Rigid Bronchoscope SET

## 5 Product Description

### 5.1 General information

#### 5.1.1 System

- Rigid Bronchoscope
- Modular system comprising a bronchoscope head (A) and replaceable tubes
- Fully dismantlable

#### 5.1.2 Bronchoscope Head



Illustration 1: Bronchoscope, components

A Bronchoscope Head

- B Bronchial / Tracheal Tube
- C Plug-in Connection
- D Injection Cannula (Luer-Lok, female)
- E Lateral Double Port
- F Axial Retainer
- G Ventilation Element
- H Axial Cap
- I Lateral Cap

- Distal end: Snap fit for connecting the tubes (B), different rotational positions possible
- Rotatable ventilation element (G) for connection to the ventilator with ISO tapered male connector (DIN EN ISO 5356-1)
- Plug-in connection (C) for lateral double port (E) or for injection cannula (D)

### 5.1.3 Tubes

*Illustration 2: Tubes*

- Replaceable tubes (B) for connection to bronchoscope head (A)
- Color coded sizes
- Proximal end: ISO tapered male connector (DIN EN ISO 5356-1) for direct connection of the respiratory tube
- Scale to measure stenosis

[ ▶ Rigid Bronchoscope Set (02BRD SET / 02BRD SET FR), page 20 ]

## 5.2 Structure and Operation

The bronchial tube mounted on the bronchoscope head is inserted through the trachea into the bronchi.

The tracheal tube mounted on the bronchoscope head is inserted into the trachea.

A telescope and instruments for diagnosis and treatment can be inserted through these tubes into the bronchi (bronchial tube) or the trachea, respectively.

To ventilate the patient the ventilation tube can be connected either to the bronchoscope head or directly to the respective tube.

### 5.3 Materials

- Stainless steel, surgical quality
- Caps: Silicon

Not made with natural rubber (latex).

No products made with natural rubber (latex) are used in the production process.

### 5.4 Other Devices to be Used in Combination with the Device

Aside from the equipment and materials required for placement, the product DUTAU-NOVATECH is intended to be used in conjunction with the following products:

Medical devices for diagnosis and therapy intended for use in rigid bronchoscopy and compatible with the length and diameter of the respective tube

## 6 Intended Use

### 6.1 Intended Purpose

The product is intended for rigid bronchoscopy.

### 6.2 Indications

Diseases / injuries that require treatment according to the intended purpose.

### 6.3 Contraindications

There are no known contraindications.

Contraindications for rigid bronchoscopy:

- Anatomical limitations of the patient, which prevent the hyperextended position of the head (e.g., cervical damage)
- Restrictions in the maxillofacial or in the laryngeal region preventing intubation using a rigid bronchoscope

### 6.4 Patient Target Group

The product is suitable for use in the following patient groups:

- Children and youth
- Adults
- Patients of all genders

### 6.5 Intended User

The intended user is a physician with experience in treating similar cases with this product or with comparable products or a physician with the following speciality:

- ENT (otorhinolaryngology)

- Pulmonology
- Thoracic surgery

The intended user will have experience with the following techniques:

- Rigid bronchoscopy

## 6.6 Expected Lifetime

[ ▶ Limitations on Processing, page 25 ]

## 6.7 Intended place of use

- Operating theatre
- Treatment room
- Endoscopic procedure room

It is the responsibility of the user to decide on a case-by-case basis which precautions must be taken for any complications that may arise.

## 7 Expected Clinical Benefit

Not applicable for this product category.

## 8 Possible Complications and Side Effects

Does not apply.

## 9 Combining with Other Procedures

Does not apply.

## 10 Shelf Life and Storage

Before first processing:

Store the product in original packaging.

## 11 Processing

The instructions listed below have been validated by the medical device manufacturer as being capable of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using equipment, materials and personnel in the processing facility achieve the desired result. This requires validation and routine monitoring of the process.

The fundamental suitability of the product for effective processing has been demonstrated by an independent, governmentally accredited and recognized test laboratory. To do so, the detergents and equipment listed in this manual were used and the procedure described in this manual was carried out.

It is possible to use detergents and equipment other than those listed in these instructions. In this case, the processor must ensure that the equipment used meets the respective criteria listed and that the processing achieves the desired result.

Observe the local regulations and hygienic rules of the medical office or hospital.

## 11.1 Warnings

### 11.1.1 General information

#### NOTICE

- Do not clean the product in the ultrasonic bath. Otherwise, the product could be damaged.
- Do not clean the product with a metal brush or steel wool. Otherwise, the product could be damaged.

For rinsing, use only low-endotoxin (maximum 0.25 endotoxin units) and sterile / low-germ (maximum 10 germs/ml) water (e.g. purified water).

### 11.1.2 Cleaning and Disinfection Agents

Do not use cleaning and disinfection agents that contain the following substances:

- Caustic soda

Corrosion inhibitors, neutralising agents and rinse aids can cause potentially critical residues on the instruments.

Do not use any rinse agents.

Only use agents suitable for cleaning / disinfection of plastic and metals.

Only use agents with proven efficacy (e.g. CE approval).

Observe all information provided by the manufacturer of the cleaning or disinfection agent (e.g. concentration, soaking time, temperature, post-rinsing).

## 11.2 Limitations on Processing

Frequent processing has little impact on these instruments. The end of the product lifetime is usually based on wear and tear as well as damage from use.

## 11.3 Cleaning Preparation

The disinfectant used for cleaning preparation is only for the safety of the user. It does not replace the disinfectant needed for cleaning and disinfection.

[ ▶ Cleaning and Disinfection Agents, page 25 ]

### 11.3.1 Preparation for Cleaning

The validation was carried out using the following equipment and method:

Cleaning solution:	Cidezyme (ASP International GmbH Zug, Switzerland), 0,5%, 20 °C
Start:	As soon as possible, but at the latest 16 h after product use
Preparation:	Open the product as widely as possible / disassemble its individual parts.
Soaking in cleaning solution:	Duration of soaking: 5 min To start the soaking time rinse all lumina at least 3 times using a single-use syringe (at least 5 ml) with cleaning solution. Brush all surfaces using a suitable brush.
Rinse:	Brush all lumina using a bottle brush in flowing water.

### 11.4 Cleaning and Disinfection

#### 11.4.1 Automated Cleaning and Disinfection

The validation was carried out using the following equipment and method:

Cleaning agent:	Neodisher Mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg), 0,5%
WD:	G 7836 CD (Miele & Cie. GmbH & Co., Gütersloh) Attach all lumina of the instrument by means of a suitable rinse adaptor to the rinsing connection of the WD.
Programme sequence:	1 min pre-cleaning in cold flowing water 5 min cleaning in 50 °C with cleaning agent solution 1 min rinse Thermal disinfection (at least 90 °C at least 5 min) Drying

#### 11.4.2 Manual Cleaning and Disinfection

Not validated for this product.

### 11.5 Drying

Use only filtered air for drying (free of oil, low in micro-organisms and particles).

### 11.6 Control, Functional Test and Maintenance

1. Inspect the product for any remaining contaminants. Clean and disinfect contaminated products again.
2. Inspect the products for damage (e.g., corrosion, damaged surfaces, deformations, other mechanical damage). Separate damaged products.
3. Maintain the threads and sliding surfaces with instrument oil (instrument oil suitable for sterilization procedure: silicone-free, based on paraffin or white oil).
4. Perform a functional check.

ATTENTION: Do not reuse discarded products.

### 11.7 Packaging

When selecting the sterilization packaging, ensure that the following requirements are met:

- Adequate protection of the product and the sterilization packaging from mechanical damage
- Adequate packaging size (tension-free sealing)
- Sterilization containers: Regular maintenance according to the manufacturer's specifications
- Suitable for steam sterilization (temperature resistant up to at least 137 °C , adequate steam permeability)

Only use standardized and certified packaging systems (EN 868 Part 2- 10, EN ISO 11607 Part 1 + 2, DIN 58953).

### 11.8 Sterilization

#### NOTICE

- The temperature in the autoclave must not exceed 137 °C. Otherwise, non-metal components of the product could be damaged.

The validation was carried out using the following equipment and method:

Sterilization procedure:	Fractionated prevacuum method
Steam sterilizer:	Selectomat HP 666-1 HRED (MMM)
Prevacuum phases:	3

Minimum temperature:	132 °C
Minimum holding time:	4 min
Minimum drying time:	20 min The actual drying time depends on parameters such as loading or sterilization parameters.

### 11.9 Storage

After sterilization, store the product dry and dust-free in the sterilization packaging.

### 11.10 Additional Information

Does not apply.

## 12 Application Instructions

### **⚠ WARNING**

- Only use intact products from intact sterile packaging or sterilization containers. Only use products that have not expired. This is the only way to ensure the product is germ-free and functional.
- Check the product and all products used in combination before and after each use for integrity, completeness and proper functioning (e.g. unintended coarse or sharp-edged surfaces, loose or incorrectly fitting parts and other obvious damage). Only use impeccable products.  
Otherwise there are risks for the health of your patient.
- Always correctly assemble the product and all products used in combination with it, including their components.  
Otherwise there are risks to the health of your patient.
- Do not bend the instruments by applying excessive pressure. Do not bend back deformed instruments.  
Otherwise, injury may be caused to the patient and user.
- When using the product: Always close unused openings of the bronchoscope head with the matching silicone caps.  
Otherwise, the ventilation of your patient is not ensured.
- Depending on the instruments used: Use silicone caps with openings to reduce the diameter of the openings of the bronchoscope head to fit.  
Otherwise, the ventilation of your patient is not ensured.

Ensure the presence of hygienic / sterile conditions needed for the intervention.

## 12.1 Product Preparation

### **⚠ WARNING**

- The product is not sterile. Process the product before first and any further application.  
This is the only way to ensure the product is germ-free and functional. [ ▶ Processing, page 24 ]

Before using the product, perform a function test.

## 12.2 Assemble and Disassemble the Bronchoscope Head

1. Push up the ventilation element against the bronchoscope head corpus.

2. Screw the axial retainer clockwise onto the bronchoscope head corpus. The axial retainer must sit tightly on the bronchoscope head corpus.

3. Push the axial cap onto the bronchoscope head to close the bronchoscope head.


To disassemble, remove all components in reverse order from the bronchoscope head.

### 12.3 Mount and Dismount Injection Cannula / Lateral Double Port

#### **⚠ WARNING**

- Ensure that the injection cannula correctly snaps into place with the plug-in connection. Otherwise, the ventilation of your patient is not ensured.
- Ensure that the jet ventilator is securely connected to the injection cannula. Otherwise, the ventilation of your patient is not ensured.

1. Push the injection cannula / the lateral double port onto the plug-in connector of the bronchoscope head as far as it will go.

2. Rotate the injection cannula against the bronchoscope head until the injection cannula audibly and noticeably snaps into place. With the lateral double port, this step is omitted.

To disassemble, pull off the injection cannula / the lateral double port from the bronchoscope head.

### 12.4 Mount the Tube onto the Bronchoscope Head and Loosen from It


Applies to the bronchial tube and the tracheal tube.

#### **⚠ WARNING**

- Ensure that the tube correctly snaps into place with the bronchoscope head. Otherwise, the ventilation of your patient is not ensured.

1. Pull back the sliding ring (A) of the bronchoscope head and hold in retracted position.

2. Introduce the tube in the desired orientation into bronchoscope head until you hear the bronchoscope head "click". Let go of the sliding ring. If the tube is introduced in the bronchoscope head correctly, the sliding ring jumps back into initial position. If needed, slightly turn the bronchoscope head.

3. Pull and rotate the tube and the bronchoscope head carefully to check that the connection is right.  
ATTENTION: Doing so, do not grasp the sliding ring of the bronchoscope head, because this would make the connection loose.

Pull the sliding ring back for disassembly. The tube will be released.

### 13 Instructing the Patient

Does not apply.

### 14 Aftercare

Does not apply.

## 15 Disposal

### WARNING

- The product was in contact with potentially infectious substances of human origin. Clean/pack the product for disposal according to the specific contamination risk. Otherwise there is a risk of infection for the user and for third parties.

Disposal must be in accordance with national disposal regulations and pursuant to the corresponding risk class.

## 16 Warranty

The reliability of the product's material and design at the time of shipment is guaranteed. The manufacturer does not know either the diagnosis of the patient or the nature of the application and has no influence on the conditions under which the product is used. The storage conditions after delivery of the product are also beyond the manufacturer's area of responsibility.

Due to biological and individual differences, no product is 100% effective under all circumstances.


**Therefore, the manufacturer cannot guarantee a positive effect or the absence of negative effects for product application. The medical staff must use the product on the basis of their medical training and experience, and they are responsible for correct application.**

The warranty (repair or replacement) applies only if the product is used in accordance with these Instructions for Use (for instruments, particularly with regard to handling, cleaning, sterilization and maintenance); the warranty period starts on the delivery date.

If you have reason to believe that a new product is faulty, please contact the Customer Service in writing immediately and provide as detailed a description as possible of the fault, the REF (product code), and the LOT (batch code) and/or series number. All allegedly defective products must be returned to us for inspection. Instruments have to be completely cleaned and sterilized, appropriate documentation must be enclosed with the return.

If the manufacturer finds that despite all due care the product was defective at the time of delivery, he will repair the product or replace it promptly. If repair or replacement of the product is not possible, the buyer has the right to cancel the purchase or to reduce the payment, but by a maximum of the purchase price amount.

Additional claims or those not mentioned here due to defect, and other claims regardless of the legal reason, including those based on illegal acts and for compensation of immaterial damages against the manufacturer, his agents, dealers and suppliers, are excluded unless existing law is contrary to the liability exclusion, e.g. in cases of intent or gross negligence or in the event of physical injury.

All claims based on the consequences of non-compliance with the Instructions for Use, including specified indications, contraindications, warnings, instructions, application, storage and off-label use, as well as the consequences of a combination with third-party products are excluded.

Furthermore, all claims that result from the use of products that have expired, or were used despite the obvious damage to the packaging, or resterilized and/or recycled contrary to the Instructions for Use, are excluded.

No one is allowed to change the above conditions, make further warranty or liability declarations, or guarantee any properties that surpass those specified in the Instructions.

The General Terms and Conditions of the manufacturer, which can be accessed at <http://www.novatech.fr/gtc> apply in all remaining instances.

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