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TWSC-UM06 Ver. 2.0 4

"Taiwan Surgical" Disposable Suction Irrigation

User Manual

PLEASE READ THE FOLLOWING INFORMATION THOROUGHLY BEFORE USING THE PRODUCT.

IMPORTANT!

These instructions for use are to assist in using this product, but not a reference to surgical techniques. All surgical procedures using the single use suction and irrigation system must be handled only by doctors who have been professionally trained and have sufficient experience to use laparoscopic and surgical techniques. Any surgeon who undertakes a procedure must understand well the risks and possible complications involved with this surgical technique. Please make sure that the patient has been notified about the manner the operation will take place, as required by local recommendations. These instructions for use are not a reference document for surgical techniques. Please refer to existing works for all information regarding surgical operations and techniques utilizing suction and irrigation devices and the relevant complications and risks. These instructions for use contain only information concerning the suction and irrigation system with and how to use it. It could not be regarded as general instructions concerning the techniques for using suction and irrigation system made by other manufactures.

INTENDED USE

Suction irrigation is indicated for general purpose use as suction irrigation device to facilitate lavage during laparoscopic surgery and general surgery. The device functions to aid flushing blood and tissue debris from the surgical site.

CONTRAINDICATIONS

- Contraindications of endoscopic techniques which apply to this product and other devices and methods. See operator's manual of your laparoscope for absolute and relative contraindications.
- Contraindications of using of this device for hysteroscopy or for cavity distention.
- Do not press irrigation key during coagulation.
- Sufficient care and distance must be maintained during use to prevent arcing to other instruments.

WARNINGS AND PRECAUTIONS

- The device is sterilized and single use. Do not reuse and re-sterilize.
- Do not use if the package has been damaged or opened.
- All surgical procedures using the single-use suction and irrigation system must be handled only by doctors who have been professionally trained and have sufficient experience to use laparoscopic and surgical technique.
- This device is obtainable only on a physician's prescription.
- Store at room temperature, 25°C. Avoid prolonged exposure to elevated temperature.
- If it is necessary to use endoscopy devices from different manufacturers in the same surgical procedure, please verify its compatibility before surgery.
- By finishing the operation, close the irrigation tube clamp, disconnect the universal flexible connector on the aspiration, and discard the whole cannula handle-tubing assembly.
- Static electricity and magnetic field influence should be noted when using the device.
- The temperatures must under 130°F (54°C) when product in transportation condition.
- The device has to be assembled on the main vacuum system of operating room, please verify its compatibility before using.
- The verified minimum and maximum pressure tolerance of the main vacuum system are as follows:

	Min	Max
Irrigation Pressure	400 mmHg (53kPa)	740 mmHg (98kPa)
Suction Pressure	400 mmHg (53kPa)	650 mmHg (86kPa)

PRODUCT DESCRIPTION

(TSI-533, TSI-545, TSI-533A, TSI-545A, DS501-4S, DSX01-3S, DSX01-4S)

The disposable suction and irrigation system consists of 2.9 meters of double tubing, with fluid bag spike and a clamp fitted on the irrigation tube, and a universal flexible connector on the suction tube. On the other side, the cannula (diameter 5 mm/10mm, length 330mm / 450mm) is connected to the tubing by the double piston hand-piece. The irrigation circuit corresponds to the blue piston, and its tubing is fitted with a clamp and the pouch spike. The aspiration circuit corresponds to the white piston, and its tubing is fitted with the universal flexible connector.

(TSI-290)

The disposable suction and irrigation system consists of 2.9 meters of double tubing, with fluid bag spike and a clamp fitted on the irrigation tube, and a universal flexible connector on the suction tube. On the other side, the cannula is connected to the Device.

SCHEMATIC VIEW

TSI-533 / TSI-545		(A) Stainless steel tube (B) Steel tube seat (C) SUCTION key (D) IRRIGATION key (E) PVC soft tube (F) Handle (G) Suction tube connector (H) Sandwicheing (I) Needle for insert (J) PVC tube (K) Needle cap
TSI-533A / TSI-545A		(A) Stainless steel tube (B) Steel tube seat (C) LOCK key (D) SUCTION key (E) IRRIGATION key (F) PVC soft tube (G) Handle (H) Suction tube connector (I) Sandwicheing (J) Needle for insert (K) PVC tube (L) Needle cap
TSI-290		(A) Suction tube connector- with the Device (B) IRRIGATION tube connector- with the Device (C) PVC soft tube (D) Suction tube connector (E) Sandwicheing (F) Needle for insert (G) PVC tube (H) Needle cap
DS501-4S/ DSX01-3S/ DSX01-4S		(A) Stainless steel tube (B) Steel tube seat

1

2

SPECIFICATIONS

Product Series	Product No.	The function of locking to continuous suction and irrigation	Stainless steel tube diameters	Stainless steel tube length	PVC tube length
TSI Series	TSI-533	--	5mm	330mm	2900mm
	TSI-545	--	5mm	450mm	2900mm
	TSI-533A	Yes	5mm	330mm	2900mm
	TSI-545A	Yes	5mm	450mm	2900mm
	TSI-290	--	--	--	2900mm
DS Series	DS501-4S	--	5mm	450mm	--
	DSX01-3S	--	10mm	330mm	--
	DSX01-4S	--	10mm	450mm	--

OPERATION NOTE

TSI Series	DS Series
Step 1. Please open the sterile packaging according to correct practice.	Step 1. Please open the sterile packaging according to correct practice.
Step 2. Please make sure the irrigation cannula, hand-piece and about one meter of tubing remain on the drape.	Step 2. Please combine stainless steel probe and make sure the irrigation cannula, hand-piece and about one meter of tubing remain on the drape.
Step 3. Please close the clamp on the irrigation tubing remove the protective cap from the spike and pierce the pouch or bottle of solution in sterile way. Aspiration of the body liquids and/or debris resulting from cleansing is effected by suction using the main vacuum system of operating room.	
Step 4. Open the irrigation tubing clamp to allow the liquid passing through. Before using the cannula, please expel all the air from the tubing. By changing a fluid bag, please close the irrigation tube clamp and replace the bag or bottle of solution in sterile method.	

3

OPERATION GUIDE

Mode	TSI-533 / TSI-545	TSI-533A / TSI-545A
1. Free mode	Free mode 	Free mode
2. Suction	Depress the SUCTION key 	Depress the SUCTION key to be locked by LOCK key
3. Irrigation	Depress the IRRIGATION key 	Depress the IRRIGATION key to be locked by LOCK key
4. Free mode	Release all the keys to back to Free mode 	Depress the LOCK key to back to Free mode

Taiwan Surgical Corporation

3F., No. 12, Sec. 2, Sheng Yi Rd., Zhubei City, Hsinchu County 30261, Taiwan
Tel: +886 3 6588129 Fax: +886 3 6588355
<http://www.twsc.com.tw>

	Catalogue number		Batch code		Sterilized using irradiation
	Date of manufacture		Manufacturer		Caution
	Use-by-date		Do not re-use		Do not re-sterilize
	Temperature limit		Keep dry		Keep away from sunlight
	Do not use if package is damaged		DEHP Free		Consult instruction for use
	EC REP	MDSG GmbH Schiffgraben 41, 30175 Hannover, Germany			

TWSC-UM06 Ver. 2.0 4